

Nextera[™] DNA Sample Preparation Kits

Sequencing's fastest and easiest sample preparation workflow, delivering libraries in 90 minutes.

Highlights

· Fastest time to results

Go from DNA to data in less than 8 hours with $MiSeq^{TM}$ system

· Easiest to use

Prepare sequencing-ready samples in 1.5 hours with 15 minutes hands-on time

Lowest DNA input

Use just 50 ng DNA per sample, enabling use with samples in limited supply

· Highest throughput

Index up to 96 samples and use master-mixed reagents to manually process > 500 samples per week



DNA to Data in Record Time

Nextera DNA Sample Preparation kits provide the fastest and easiest workflow, enabling sequencing-ready libraries to be generated in less than 90 minutes, with less than 15 minutes of hands-on time. DNA is simultaneously fragmented and tagged with sequencing adapters in a single step, using standard lab equipment. Libraries prepared with Nextera kits are compatible with all Illumina® sequencers (Table 1).

Table 1: Nextera DNA Sample Prep Specifications

Specification	Value	
Input DNA	50 ng	
Available indices	Up to 96	
Compatible sequencers	HiSeq® 1000/2000, HiScanSQ®, Genome Analyzer IIx and MiSeq systems	
Read lengths supported	Supports all read lengths on any Illumina sequencing system, including 2 x 150 reads on the MiSeq System	
Typical median insert size	~250 bp	
Sample DNA input type	Genomic DNA and PCR amplicons	


Breakthrough Chemistry

Nextera technology employs a single "tagmentation" reaction to simultaneously fragment and tag DNA with adapters (Figure 2). This process occurs in a single step using master-mixed reagents to provide PCR-ready templates in as little as 15 minutes. Sequencing adaptors and indices are then added to the gDNA fragment by PCR. The optimized Nextera PCR protocol leads to improved performance with GC regions. From start to finish, the complete Nextera sample preparation protocol is over 80% faster than any other method available.

Improved Multiplexing

Nextera DNA Sample Preparation kits feature an innovative indexing solution for processing and uniquely barcoding up to 96 samples. Multisample studies can be conveniently managed using the Illumina Experiment Manager, a freely available software tool that provides easy reaction setup for plate-based processing.

Following the addition of two indices to each gDNA fragment, up to 96 uniquely indexed samples can be pooled and sequenced together in a single lane on any Illumina sequencer. After sequencing, the unique combination of the two indices is used to demultiplex the data and assign reads to the proper sample in the pool. Using this dual barcode approach, Nextera Index Kits only require 20 unique index oligos to process up to 96 samples, making this an easily scalable approach for sample indexing.



Accelerated Applications

Nextera DNA Sample Preparation kits are ideal for experiments where speed and ease are paramount. The low 50 ng DNA input also makes this method amenable to precious samples available in limited quantity. This sample preparation workflow can shorten the overall sequencing workflow time for a wide variety of established applications¹⁻⁷ and can be automated for even greater throughput. The combination of the MiSeq system and Nextera sample preparation kits provide rapid DNA to data in as little as 8 hours, enabling rapid applications such as small genome and amplicon sequencing (Table 2).

Summary

Illumina's Nextera DNA Sample Preparation kit provides sequencing's fastest and easiest sample preparation workflow, delivering completed libraries in 90 minutes that are compatible with all Illumina sequencing systems. Nextera enables high-throughput studies with a built-in solution for indexing up to 96 samples with ultra low DNA input. Combined with the MiSeq system, Nextera DNA Sample Preparation Kits enable the fastest DNA to data—all in a single day.

Table 2: Representative Nextera Applications

Examples of Nextera Applications		
Small-genome resequencing		
Small-genome de novo sequencing		
Amplicon sequencing		
Clone or plasmid sequencing		

References

- Ramirez MS, Adams MD, Bonomo RA, Centrón D, et al. (2011) Genomic analysis of *Acinetobacter baumannii* A118 by comparison of optical maps: Identification of structures related to its susceptibility phenotype. Antimicrob Agents Chemother, 55(4): 1520–6.
- Adey A, Morrison HG, Asan, Xun X, Kitzman JO, et al. (2010) Rapid, low-input, low-bias construction of shotgun fragment libraries by high-density in vitro transposition. Genome Biol 11: R119.
- Bimber BN, Dudley DM, Lauck M, Becker EA, Chin EN, et al. (2010)
 Whole-genome characterization of human and simian immunodeficiency virus intrahost diversity by ultradeep pyrosequencing. J Virol 84: 12087–92.
- Kitzman JO, Mackenzie AP, Adey A, Hiatt JB, Patwardhan RP, et al. (2010) Haplotype-resolved genome sequencing of a Gujarati Indian individual. Nat Biotechnol 29: 59–63.
- Linnarsson, S. (2010) Recent advances in DNA sequencing methods -General principles of sample preparation. Exp Cell Res 316: 1339–43.
- Sudmant PH, Kitzman JO, Antonacci F, Alkan C, Malig M, et al. (2010) Diversity of human copy number variation and multicopy genes. Science 330: 641–646.
- Voelkerding KV, Dames S, and JD Durtschi (2010) Next generation sequencing for clinical diagnostics-Principles and application to targeted resequencing for hypertrophic cardiomyopathy. J Mol Diagn 12: 539–551.

Ordering Information

Product	Catalog No.	
Nextera DNA Sample Preparation Kit (96 Samples)	FC-121-1031	
Nextera DNA Sample Preparation Kit (24 Samples)	FC-121-1030	
Nextera Index Kit (96 Indices, 384 Samples)	FC-121-1012	
Nextera Index Kit (24 Indices, 96 Samples)	FC-121-1011	
TruSeq Dual Index Sequencing Primer Kit, Single Read (single-use kit)	FC-121-1003	
TruSeq Dual Index Sequencing Primer Kit, Paired-End Read (single-use kit)	PE-121-1003	

Illumina, Inc. • 9885 Towne Centre Drive, San Diego, CA 92121 USA • 1.800.809.4566 toll-free • 1.858.202.4566 tel • techsupport@illumina.com • illumina.com

FOR RESEARCH USE ONLY

© 2011 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Pub. No. 770-2011-021 Current as of 07 October 2011

